

## Solvent preparation



### Step 1

- Get fresh THF from the solvent system
- If your THF is already degassed, you can skip this step



### Step 2

- You will need a high purity nitrogen tank or argon tank
- Open the schlenk flask using standard schlenk techniques
- Submerge the long needle in the THF solution, and add a vent needle (you should see bubbles now)
- Sonicate for 30 minutes
- THF is now degassed

## Schlenk preparation



### Step 3

- Oven dry a 100 mL Schlenk overnight with a stir bar
- You will need two spatulas, and a funnel that has a long neck




### Step 4

- Find Sm metal and diiodoethane
- You need 1.65 g of Sm
- You need 1.55 g of diiodoethane
- This makes 55 mL of 0.2 M solution of  $\text{SmI}_2$



## Making the 0.2 M SmI<sub>2</sub> solution



### Step 5

- Bring the filled Schlenk back to your fume hood
- Connect to vacuum and cycle three times



### Step 6

- Add 55 mL of degassed THF
- Turn stirring on



### Step 7

- Ethylene is generated (bubbles circled in red), so pull vacuum on the Schlenk, but turn the stopcock quickly otherwise SmI<sub>2</sub> will be pulled up into the Schlenk line
- The solution will be warm to the touch



### Step 8

- Stir for 20 minutes while open to the nitrogen line, the solution should stop bubbling, and cool to room temperature
- Now the Schlenk can be sealed and left to stir overnight
- Store in freezer
- Enjoy your 0.2 M SmI<sub>2</sub>!

